

Cytology
**endocervical
endometrial
brushes**

Distributed by



Eurogine, S.L.
C/Antonio Gaudí, 91
08830 Sant Boi de Llobregat, Barcelona (Spain)
Tel.: +34 93 630 43 45 - Fax: +34 93 630 39 71

Web: www.eurogine.com
international@eurogine.com

Brushes to collect endocervical / endometrial cells

Ref.	Description	Packaging
07010000	Endocervical brush (bags of 10 units)	Box of 100 units
07010100	Endocervical brush (bags of 10 units)	Box of 500 units
07020000	Sterile endocervical brush	Box of 100 units
07030000	Sterile endometrial brush	Box of 100 units



Miscellaneous

Ref. 11500000 PUNCTURE DEVICE FOR BIOPSY

This device designed by Eurogine permits the use of disposable 10cc and 20cc syringes providing almost their total capacity of vacuum aspiration. It allows having one hand free. Ambidextrous design.

Packaging / Unit



Miscellaneous

Ref. 18000000 NON-SCALPEL VASECTOMY KIT

Kit containing two different Ø ring fixation clamps and one sharpened mosquito forceps. Stainless steel. NSV has several advantages with the same effectiveness as traditional vasectomy procedure:

- Minimizes pain, infections and bleeding
- No stitching is needed
- Reduces time of surgical intervention

Packaging / Unit



Ref. 11260200 VAGINAL DILATORS FOR VAGINISMUS

Case containing six therapeutic vaginal dilators of gradually increasing sizes and a lubricant for the treatment of vaginismus.

Packaging / Unit



Ref. 1857-290 STAINLESS STEEL FORCEPS FOR IUDS REMOVING

Stainless steel forceps developed by Eurogine for IUD removing or arranging retracted threads.

Packaging / Unit



Eurogine, S.L.

C/Antonio Gaudí, 91
08830 Sant Boi de Llobregat, Barcelona (Spain)
Tel.: +34 93 630 43 45 - Fax: +34 93 630 39 71

Web: www.eurogine.com
international@eurogine.com

Intrauterine devices for contraception

The largest known range of IUD's in the international market

Gold T 375

COPPER+GOLD (Cu375+Au)

- 99,9% pure copper wire + gold core Ø 0,1mm
- Y-shaped (NOVAPLUS)
- 3 sizes available
- Unique and easy insertion system

Ref.	Model	Size
01040000	Gold T 375	Maxi (Large)
01040100	Gold T 375	Normal (Medium)
01040200	Gold T 375	Mini (Small)



COPPER + SILVER (Cu+Ag)

- 99,9% pure copper wire + silver core Ø 0,1mm
- Y-shaped (NOVAPLUS) or α -shaped depending on the model
- 1, 2 or 3 sizes available depending on the model
- Cu200+Ag, Cu250+Ag, Cu375+Ag, Cu380+Ag

Ref.	Model	Size
01010500	Novaplug T 380 Ag	Normal (Medium)
01010600	Novaplug T 380 Ag	Mini (Small)
01010700	Novaplug T 380 Ag	Maxi (Large)
01010000	T 200 Ag	Normal (Medium)
01010100	T 200 Ag	Mini (Small)
01030300	Ancora 250 Ag	Normal (Medium)
01030400	Ancora 375 Ag	Normal (Medium)

COPPER (Cu)

- 99,9% pure copper wire
- Y-shaped (NOVAPLUS) or α -shaped depending on the model
- 1, 2 or 3 sizes available depending on the model
- Cu200, Cu250, Cu375, Cu380

Ref.	Model	Size
01020000	T 200 Cu	Normal (Medium)
01030100	Ancora 250 Cu	Normal (Medium)
01030200	Ancora 250 Cu	Mini (Small)
01030000	Ancora 375 Cu	Normal (Medium)
01020100	Novaplug T 380 Cu	Normal (Medium)
01020200	Novaplug T 380 Cu	Mini (Small)
01050000	Copper T 380A	Normal (Medium)



Ref. 03020000 BIASPIR

Sterile disposable device for endometrial sampling via double aspiration.

- Cannula with 2mm holes for aspiration
- Cannula with 3mm holes for aspiration
- Vacuum syringe with lock system
- 12ml plastic tube for sample collection
- Patient Id. label

Packaging / Box of 25 units

Disposable kits for endometrial sampling



Ref. 03010000 GINEASPIR

Sterile disposable device for endometrial sampling via curettage and vacuum.

- 4mm flexible cannula
- Vacuum syringe with lock system
- 12ml plastic tube for sample collection
- Patient Id. label

Packaging / Box of 25 units

Ref. 03030000 CYTOSPAT

Sterile disposable device for endometrial sampling via curettage.

- Plastic curette with a rhombus end / blade effect
- Vacuum syringe with lock system
- 12ml plastic tube for sample collection
- Patient Id. label

Packaging / Box of 25 units

